










APC (All Purpose Cuff) Infant / Child / Small Adult / Adult / Large Adult / Thigh Size:

	Infant Range 8 – 13 cm
	Child Range 12-19cm Length 30cm
	Small Adult Range 17-25cm Length 43cm
	Small Adult Long Range 17-25cm Length 50cm
	Adult Range 23-33cm Length 50cm
	Adult Long Range 23 - 33 cm Length 63 cm



	Large Adult Range 31 - 40 cm Length 64 cm
	Large Adult Long Range 31 - 40 cm Length 80 cm
	Thigh Range 38 - 50 cm